

Raw sequences files
and quality control

Fastq Format

This format is designed to handle base quality metrics output from sequencing machines.

```
Identifier  ———| @HWI-EAS209_0006_FC706VJ:5:58:5894:21141#ATCACG/1
Sequence   ———| TTAATTGGTAAATAAATCTCCTAATAGCTTAGATNTTACCTTNNNNNNNNNNNTAGTTTCTTGAGA
+ sign & identifier ———| +HWI-EAS209_0006_FC706VJ:5:58:5894:21141#ATCACG/1
Quality scores ———| efcfffffcfeeffffcfffffdddf`feed)`]_Ba_^__[YBBBBBBBBBBBRTT\]][] dddd`
```

Base T
phred Quality] = 29

Line 1 begins with the '@' character and is followed by a sequence identifier and an optional description.

Line 2 is the sequence letters.

Line 3 begins with a '+' character; it marks the end of the sequence and is optionally followed by the same sequence identifier again in line 1.

Line 4 encodes the quality values for the sequence in Line 2, and must contain the same number of symbols as letters in the sequence.

Quality scores

@EAS139:136:FC706VJ:2:2104:15343:197393 1:Y:18:ATCACG

CCGTCAATTCATTAGTTTTAACCCTTGCGGCCGTACTCCCCAGGCGGT

+

AAAAAAAAAAAAA:9@: : : ? ? @ @ : : FF AAAA ACCAA: : : : BB @ @ ? A ?

ASCII encoding

40:@	90:Z	141:a
41:A	91:[142:b
42:B	92:\	143:c
43:C	93:]	144:d
44:D	94:^	145:e
45:E	95:_	146:f
... :...	... :...	... :...


$$\text{Phred} = -10 \log_{10} p$$

p = Probability call is incorrect


Quality Score	Probability of incorrect base call	Base call accuracy
10	1 in 10	90%
20	1 in 100	99%
30	1 in 1000	99.9%
40	1 in 10000	99.99%

Assess quality with FastQC


Per base sequence quality




Per tile sequence quality




Per sequence quality scores




Per base sequence content



Per sequence GC



Sequence length



<https://training.galaxyproject.org/training-material/topics/sequence-analysis/tutorials/quality-control/tutorial.html#sequence-length-distribution>

Let's have a look at the first
few sequences and check the
sequencing quality with fastqc